

My Bike® Program Order Form:

Please note that this form should only be completed once the child's application is approved by Variety, and it should only be done by the child's physical / occupational therapist, primary care physician, physician assistant, or certified nurse practitioner.

The person completing the form must do so while the child is on a Rifton adaptive bike.

If a Rifton adaptive bike is not available in your community, then a My Bike® fitting/ evaluation will need scheduled by Variety.



Child's Name:	Date:
Name of Person Completing Form:	Title/Affiliation:
Phone Number and Email Address:	
Signature of Person Completing Form:	

	Name:	Specifics:	Small (R120):	Medium (R130):	Large (R140):
	Tricycle Frame includes self-leveling pedals with straps, backrest tube with handle, tote box for storage, and handbrake for the large only.		R112 □	R113 □	R114 □
Required Components	Color (Choose one color)		Blue ☐ Red ☐ Yellow ☐	Blue ☐ Red ☐ Yellow ☐	Blue ☐ Red ☐ Yellow ☐
			Pink Lime	Pink Lime	Pink ☐ Lime ☐
	Seat gg	Small	R122 □	R132 □	N/A
	Seat (with seatbelt)	Large	R128 □	R138 □	(included)
	Handlebar 요.	Straight	R126 □	R136 □	R146 🗖
		Loop	R129 □	R139 □	R149 🗖
	Back Support	Trunk Support System	R147 □	R147 □	R147 □
	Back Support SO	Backrest Pad	R145 □	R145 🗖	R145 🗖
	Headrest	Winged	R159 🗖	R159 □	R159 🗖
	(requires trunk support system)	Contoured	R156 □	R156 □	R156 □
Optional Accessories	ř	Flat	R157 🗖	R157 □	R157 🗖
	Handbrake (Must be factory installed. Cannot be added later)		Not Available	R125 □	(included)
	Communication Tray	Small	K225 🗖	K225 □	K225 🗖
	Communication Tray	Large	N/A	K245 □	K245 □
	Abductor		R121 □	R131 □	R141 □
	Front Guide Bar		R135 🗖	R135 □	R135 🗖
	Rear Steering Bar		R123 □	R123 🗖	R123 🗖
	Front Pulley		R144 □	R144 □	R144 □
	Additional Backrest Pad		R145 □	R145 □	R145 🗖
	Stationary Stand		(included by Variety)	(included by Variety)	(included by Variety)



If you need help or have questions, please call Variety at 724-933-0460.



The extended leg should reach from seat to pedal comfortably when both feet and torso are strapped into place.

Important: User's weight must not exceed the maximum working load. Rifton tricycles are not intended for clients with strong reflexes or poorly developed trunk balance and head control. Consult your client's therapist before ordering.

Item dimensions (inches)	R120 small	R130 medium	R140 large
Overall height (at lowest position)	281/2	31	361/2
Overall length	46	52	65
Overall width	271/2	27½	31½
Wheel diameter	13	15	19
Seat center to extended pedal	17-22	22 – 28	25 – 35
Seat to top of trunk support	15½ – 22	15½-22	15½ –22
Gear ratio (approximate)	1:1	1:1	1:1
Item weight (lbs)	47	56	65
Distance between trunk laterals	8-14	8 – 14	8 – 14
Max. working load (lbs)	150	160	200